

GENIUSlink Basic Information

Charging Time

- 2-hour maximum (no power to full power)
- Charge nightly for best performance

Wearing Time

- Continuous Bluetooth streaming (phone or audio): 4-6 hours
- Continuous audio streaming via mini jack port: up to 10 hours
- Remote Control use only: up to 1 week
- For wearing time beyond limits GENIUSlink can be used while plugged into an AC source or to a portable charger.

Preventative Maintenance

- Using the silicone sleeve helps to prevent moisture and dirt from entering open ports/connections of the GENIUSlink Medallion.

Setup Sequence

1. Charge the GENIUSlink
2. Perform a first fit for hearing aids
3. Pair hearing aids to GENIUSlink
4. Pair GENIUSlink to phone

LED Light Patterns

Flashing lights indicate GENIUSlink is **searching** for something

- **Green flashing:** searching for hearing aids
 - **Blue flashing:** searching to pair to a Bluetooth device (phone, tablet, etc.)
 - **Green and Blue flashing simultaneously:** first time pairing-searching for both hearing aids and Bluetooth devices.
 - **Green and Blue flashing alternately:** searching to pair to a Transmitter/MEComConnect
- *Note: Orange flashing:** error - connection lost with hearing aids (only visible when MFB or volume up/down pressed), neckloop connection loose, or needs recharging

Solid lights indicate Smart Connect is actively **connected** and/or streaming

- **Blue solid:** transmitting phone call or streaming audio to hearing aids.
- **Green solid:** only visible for 2-3 seconds when first connecting to hearing aids, or when plugged in and fully charged
- **Orange solid:** only visible when plugged in and charging.