Use myMEMSI for fast and accurate service.

///\/\/\/Miracle-Ear	[®] July 2017		INIGED GROEK I GRAN
	MM / DD .	/ YY TEL	: (800) 314-2694 FAX: (800) 234-0980
Account #	Date: / /	P.O. #	
Ship To Address:	Contact:		
	Tel:		
Client Information Please use capital letter	s and print neatly.	Hearing Loss	
Last Name, First Name, Middle Initial		L R	L R
		Mild	□ □ Severe
		☐ ☐ Moderate	□ □ Profound
Shipping 2nd Day (Standard)	Overnight Shipping		
Earmolds will ship with the BTE when a			
STEP 1 - Material		TEP 5 - Tubing Size / Type	STEP 7 - Removal String
L R	_	13 Tube	L R
Acrylic (Hard)	Gloss Matte L R L R	Regular or □ Double Bend L R	Add Removal String* *Standard on all RIC Earmolds,
Silicone (Soft)	Clear	3.1mm OD	Open, and Open with Lock
STEP 2 - Style ¹ Hard material only		3.3mm OD	Repair & Remake
Receiver-in-the-Canal (RIC) Earmolds choose receiver size and style L R L R L R		Ory Tube 3.3mm OD	Earmold Serial Number
GENIUS™ 2.0, S M P		Quick-Connect 3.1mm OD	
GENIUS™/ClearVation (X400, X300, X200 & X175) □ □ □ □ □ □	material with gloss imish	Quick-Connect Dry 3.3mm OD	
ClearAmp 45 dB 55 dB 65 dB	Red Opaque 🔲 🔲	Quick-Connect 3.6mm OD Tube Lock is the default connection	R
(X150)	Orange Opaque 🔲 🔲	for all silicone earmolds	Repair - All out of warranty repairs
RIC (AR)	Tellow Opaque L L L	Open Tube*	will incur a charge The earmold will be evaluated
Foil RIC (FR) ¹	Blue Opaque □ □	1 2 3 4 5 6	and repaired if possible.
Foil RIC with Lock (FL) ¹	Lilac Opaque	ielect a Length L 🔲 🗎 🗎 🗎 🗎	Replace receiver adapter
Canal RIC (AG)	Black Opaque	Open Tube (for X300 PX, X200 M,	Add/Replace elbow/anchor
Canal with Lock RIC (AH)		(200 P, X200 PX)*	Increase/Decrease tube size
Half-Shell RIC (AB)		1 2 3 4 5 6 Select a Length L 🗆 🗆 🗆 🗆	Add/Replace removal option
Skeleton RIC (AC) 0	Blue Translucent 🔲 🔲	R	Increase/Decrease vent size
1/2 Skeleton RIC (AF) 1 \square \square \square \square \square \square \square Semi-Skeleton RIC (AD) 1 \square \square \square \square \square \square \square \square	Purple Translucent	(For GENIUS™ 2.0/GENIUS SX only) Thin Tube* (For GENIUS™ 2.0 SX Only)	Remake*
Full Shell RIC (AA)	Red Glitter	1 2 3 4 5** Select a Length L	90-day remake allowance from original invoice date. *Charges may apply.
Skeleton Open RIC (AS) ¹	Yellow Glitter	R 🗆 🗆 🗆	Reason for remake:
Standard Earmolds L R Canal (CT)	Green Glitter	*Available with any Acrylic earmold **Length 5 must be glued	L R
Canal with Lock (CL)	Blue Glitter	TEP 6 - Venting	Too loose
Skeleton (SK) \square \square Skeleton Open (SO) ¹ \square \square	Purple Glitter	L R	Feedback 🗆 🗆
Semi-Skeleton (SS) ¹	Clear Clitter	argest Vent Possible	Canal too short
3/4 Skeleton (FC)¹ □ □ 1/2 Skeleton (HC)¹ □ □	Swirl*	Semi-IROS Gelect-a-Vent 3.0mm	Manufacturer error
Half-Shell (HS)	choice in the special instructions section	Select-a-Vent 2.0mm	Earmold sticks out
Full-Shell (SH)	Defaulte in Peld	Standard 1.6mm	Wire/tube protrudes
Open Earmolds		Power Vent 1.2mm*	In-Office Use
Open (LI) ¹ □ □ Open with Lock (LS) ¹ □ □	N	No Vent	L R Use modified earmold/ □ □
Choose an OpenTube length. Vent Selection N/A.		Skeleton and Full-Shell earmolds	impression for the remake
STEP 4 - Canal Length			Special Instructions:
L R Deep	Deep canal Impression should be 2mm		
Long	Medium canal past the 2nd bend to provide	5mm	
Short	the best cosmetic wire/	To build any receiver-in-the- canal earmold the canal	
Cut to Mark	tube placement	must be at least 5mm wide at the narrowest point	

©2017 Miracle-Ear, Inc. 14452FORM/F